

### **EDTA - concentrate to make 1L 0.1M Solution**

Art. ID	LM24-OP-2001955
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

#### **Description**

Concentrated Volumetric Solutions / Volumetric concentrates offer the convenience of prepared volumetric solutions with the benefits of reduced shipping costs and storage space. Concentrated are designed to be quantitatively transferred and brought to volume. / Certificate of Analysis with actual contents, lot number, expiry date and traceability to SI of NIST (SRM) or of Laboratories/ Producers accredited according to ISO/ IEC 17025 and/ or ISO Guide 34 / Concentrated volumetric solutions in sealed bottles. Easy to use. No ampoules to break. No need for pipetting. Empty the bottle in a volumetric flask, rinse and complete to volume /