

Potassium Permanganate concentrate to make 1L 1/50M (0.1N) KMnO4 (DRUG PRECURSOR 2)

Art. ID LM24-OP-2001970
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

Concentrated Volumetric Solutions / Volumetric concentrates offer the convenience of prepared volumetric solutions with the benefits of reduced shipping costs and storage space. Concentrated are designed to be quantitatively transferred and brought to volume. / Certificate of Analysis with actual contents, lot number, expiry date and traceability to SI of NIST (SRM) or of Laboratories/ Producers accredited according to ISO/ IEC 17025 and/ or ISO Guide 34 / Concentrated volumetric solutions in sealed bottles. Easy to use. No ampoules to break. No need for pipetting. Empty the bottle in a volumetric flask, rinse and complete to volume /

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Contains a DRUG PRECURSOR 2(B)	Potassium Permanganate	[7722-64-7]				