

Sodium Hydroxide - concentrate to make 1L 1M (1N) Solution

| | |
|-----------------|--|
| Art. ID | LM24-OP-2001972 |
| Unit | 100 mL |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1824 II E2 |

Description

Concentrated Volumetric Solutions / Volumetric concentrates offer the convenience of prepared volumetric solutions with the benefits of reduced shipping costs and storage space. Concentrated are designed to be quantitatively transferred and brought to volume. / Certificate of Analysis with actual contents, lot number, expiry date and traceability to SI of NIST (SRM) or of Laboratories/ Producers accredited according to ISO/ IEC 17025 and/ or ISO Guide 34 / Concentrated volumetric solutions in sealed bottles. Easy to use. No ampoules to break. No need for pipetting. Empty the bottle in a volumetric flask, rinse and complete to volume /