

### **Sodium Hydroxide -concentrate to make 1L 0.1M (0.1N) Solution**

Art. ID	LM24-OP-2001973
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1824 III E1

#### **Description**

Concentrated Volumetric Solutions / Volumetric concentrates offer the convenience of prepared volumetric solutions with the benefits of reduced shipping costs and storage space. Concentrated are designed to be quantitatively transferred and brought to volume. / Certificate of Analysis with actual contents, lot number, expiry date and traceability to SI of NIST (SRM) or of Laboratories/ Producers accredited according to ISO/ IEC 17025 and/ or ISO Guide 34 / Concentrated volumetric solutions in sealed bottles. Easy to use. No ampoules to break. No need for pipetting. Empty the bottle in a volumetric flask, rinse and complete to volume /