

Caffeine 0.15 mg/mL in Water (ISO 17034 \ ISO 17025)

Art. ID LM24-OP-2002060
Unit 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) /not restricted

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | Caffeine | [58-08-2] | | | | |