

Supelco 37 Component FAME Mix, certified reference material, TraceCERT®, in dichloromethane (varied conc.), ampule of 1 mL

Art. ID SAF-CRM47885
Unit 1 mL
Delivery details Dry ice shipment UN1593 III E1

Description

TraceCERT® is a trademark of Sigma-Aldrich

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methyl butyrate	[623-42-7]	400	µg/mL		
	Methyl hexanoate	[106-70-7]	400	µg/mL		
	Methyl octanoate	[111-11-5]	400	µg/mL		
	Methyl decanoate	[110-42-9]	400	µg/mL		
	Methyl undecanoate	[1731-86-8]	200	µg/mL		
	Methyl laurate (C12:0)	[111-82-0]	400	µg/mL		
	Methyl tridecanoate	[1731-88-0]	200	µg/mL		
	Methyl myristate (C14:0)	[124-10-7]	400	µg/mL		
Methyl myristoleate	Myristoleic acid methyl ester	[56219-06-8]	200	µg/mL		
	Methyl Pentadecanoate	[7132-64-1]	200	µg/mL		
	cis-10-Pentadecenoic acid methyl ester	[90176-52-6]	200	µg/mL		
	Methyl palmitate	[112-39-0]	600	µg/mL		
	Methyl palmitoleate	[1120-25-8]	200	µg/mL		
	Methyl heptadecanoate	[1731-92-6]	200	µg/mL		
	Methyl cis-10-heptadecenoate	[75190-82-8]	200	µg/mL		
	Methyl stearate	[112-61-8]	400	µg/mL		
	9-trans-Octadecenoic acid methyl ester (Elaidic acid methyl ester)	[1937-62-8]	200	µg/mL		
cis-9-Oleic acid methyl ester	(Z)-9-Octadecenoic acid methyl ester (Oleic acid methyl ester)	[112-62-9]	400	µg/mL		
Methyl linolelaidate	Methyl linoelaidate (9,12tr), No. of C: C18:2	[2566-97-4]	200	µg/mL		

	Methyl linoleate (C18:2)	[112-63-0]	200	µg/mL
	Methyl arachidate (C20:0)	[1120-28-1]	400	µg/mL
Methyl gamma-linolenate	6,9,12-Octadecatrienoic acid methyl ester (gamma-Linolenic acid methyl ester)	[16326-32-2]	200	µg/mL
	Methyl cis-11-eicosenoate	[2390-09-2]	200	µg/mL
	Methyl linolenate (C18:3)	[301-00-8]	200	µg/mL
	Methyl heneicosanoate	[6064-90-0]	200	µg/mL
cis-11,14-Eicosadienoic acid methyl ester	Methyl eicosadienoate	[2463-02-7]	200	µg/mL
	Methyl behenate (C22:0)	[929-77-1]	400	µg/mL
cis-8,11,14-Eicosatrienoic acid methyl ester	Methyl cis-8,11,14-eicosatrienoate	[21061-10-9]	200	µg/mL
	Methyl erucate (C22:1)	[1120-34-9]	200	µg/mL
	Methyl 11-14-17-eicosatrienoate (C20:3)	[55682-88-7]	200	µg/mL
	Methyl tricosanoate (C23:0)	[2433-97-8]	200	µg/mL
	5,8,11,14-Eicosatetraenoic acid methyl ester (Arachidonic acid methyl ester)	[2566-89-4]	200	µg/mL
	Methyl cis-13,16-docosadienoate	[61012-47-3]	200	µg/mL
	Methyl lignocerate (C24:0)	[2442-49-1]	400	µg/mL
	cis-5,8,11,14,17-Eicosapentaenoic acid methyl ester	[2734-47-6]	200	µg/mL
	Methyl nervonate (C24:1)	[2733-88-2]	200	µg/mL
	Methyl docosahexaenoate (C22:6)	[2566-90-7]	200	µg/mL

