

CONOSTAN - Multi-element standard in oil S-12+Mo: 200 ppm

Art. ID CON-150-012-021
Unit 200 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,02 | % | | |
| | Aluminium (Al) | [7429-90-5] | 200 | µg/g | | |
| | Chromium (Cr) | [7440-47-3] | 0,02 | % | | |
| | Chromium (Cr) | [7440-47-3] | 200 | µg/g | | |
| | Copper (Cu) | [7440-50-8] | 0,02 | % | | |
| | Copper (Cu) | [7440-50-8] | 200 | µg/g | | |
| | Iron (Fe) | [7439-89-6] | 0,02 | % | | |
| | Iron (Fe) | [7439-89-6] | 200 | µg/g | | |
| | Lead (Pb) | [7439-92-1] | 0,02 | % | | |
| | Lead (Pb) | [7439-92-1] | 200 | µg/g | | |
| | Magnesium (Mg) | [7439-95-4] | 0,02 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 200 | µg/g | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,02 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 200 | µg/g | | |
| | Nickel (Ni) | [7440-02-0] | 0,02 | % | | |
| | Nickel (Ni) | [7440-02-0] | 200 | µg/g | | |
| | Silicon (Si) | [7440-21-3] | 0,02 | % | | |
| | Silicon (Si) | [7440-21-3] | 200 | µg/g | | |
| | Silver (Ag) | [7440-22-4] | 0,02 | % | | |
| | Silver (Ag) | [7440-22-4] | 200 | µg/g | | |
| | Sodium (Na) | [7440-23-5] | 0,02 | % | | |
| | Sodium (Na) | [7440-23-5] | 200 | µg/g | | |
| | Tin (Sn) | [7440-31-5] | 0,02 | % | | |
| | Tin (Sn) | [7440-31-5] | 200 | µg/g | | |
| | Titanium (Ti) | [7440-32-6] | 0,02 | % | | |
| | Titanium (Ti) | [7440-32-6] | 200 | µg/g | | |