

## **CONOSTAN - Multi-element standard in oil S-12+Mo: 500 ppm**

Art. ID                      CON-150-012-027  
Unit                         200 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,05	%		
	Aluminium (Al)	[7429-90-5]	500	µg/g		
	Chromium (Cr)	[7440-47-3]	0,05	%		
	Chromium (Cr)	[7440-47-3]	500	µg/g		
	Copper (Cu)	[7440-50-8]	0,05	%		
	Copper (Cu)	[7440-50-8]	500	µg/g		
	Iron (Fe)	[7439-89-6]	0,05	%		
	Iron (Fe)	[7439-89-6]	500	µg/g		
	Lead (Pb)	[7439-92-1]	0,05	%		
	Lead (Pb)	[7439-92-1]	500	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,05	%		
	Magnesium (Mg)	[7439-95-4]	500	µg/g		
	Molybdenum (Mo)	[7439-98-7]	0,05	%		
	Molybdenum (Mo)	[7439-98-7]	500	µg/g		
	Nickel (Ni)	[7440-02-0]	0,05	%		
	Nickel (Ni)	[7440-02-0]	500	µg/g		
	Silicon (Si)	[7440-21-3]	0,05	%		
	Silicon (Si)	[7440-21-3]	500	µg/g		
	Silver (Ag)	[7440-22-4]	0,05	%		
	Silver (Ag)	[7440-22-4]	500	µg/g		
	Sodium (Na)	[7440-23-5]	0,05	%		
	Sodium (Na)	[7440-23-5]	500	µg/g		
	Tin (Sn)	[7440-31-5]	0,05	%		
	Tin (Sn)	[7440-31-5]	500	µg/g		
	Titanium (Ti)	[7440-32-6]	0,05	%		
	Titanium (Ti)	[7440-32-6]	500	µg/g		