

CONOSTAN - Multi-element standard in oil S-21+K+Li: 500 ppm

Art. ID CON-150-021-064
Unit 200 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,05	%		
	Aluminium (Al)	[7429-90-5]	500	µg/g		
	Barium (Ba)	[7440-39-3]	0,05	%		
	Barium (Ba)	[7440-39-3]	500	µg/g		
	Boron (B)	[7440-42-8]	0,05	%		
	Boron (B)	[7440-42-8]	500	µg/g		
	Cadmium (Cd)	[7440-43-9]	0,05	%		
	Cadmium (Cd)	[7440-43-9]	500	µg/g		
	Calcium (Ca)	[7440-70-2]	0,05	%		
	Calcium (Ca)	[7440-70-2]	500	µg/g		
	Chromium (Cr)	[7440-47-3]	0,05	%		
	Chromium (Cr)	[7440-47-3]	500	µg/g		
	Copper (Cu)	[7440-50-8]	0,05	%		
	Copper (Cu)	[7440-50-8]	500	µg/g		
	Iron (Fe)	[7439-89-6]	0,05	%		
	Iron (Fe)	[7439-89-6]	500	µg/g		
	Lead (Pb)	[7439-92-1]	0,05	%		
	Lead (Pb)	[7439-92-1]	500	µg/g		
	Lithium (Li)	[7439-93-2]	0,05	%		
	Lithium (Li)	[7439-93-2]	500	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,05	%		
	Magnesium (Mg)	[7439-95-4]	500	µg/g		
	Manganese (Mn)	[7439-96-5]	0,05	%		
	Manganese (Mn)	[7439-96-5]	500	µg/g		
	Molybdenum (Mo)	[7439-98-7]	0,05	%		
	Molybdenum (Mo)	[7439-98-7]	500	µg/g		
	Nickel (Ni)	[7440-02-0]	0,05	%		
	Nickel (Ni)	[7440-02-0]	500	µg/g		
	Phosphorus (P)	[7723-14-0]	0,05	%		
	Phosphorus (P)	[7723-14-0]	500	µg/g		
	Potassium (K)	[7440-09-7]	0,05	%		

Potassium (K)	[7440-09-7]	500	µg/g
Silicon (Si)	[7440-21-3]	0,05	%
Silicon (Si)	[7440-21-3]	500	µg/g
Silver (Ag)	[7440-22-4]	0,05	%
Silver (Ag)	[7440-22-4]	500	µg/g
Sodium (Na)	[7440-23-5]	0,05	%
Sodium (Na)	[7440-23-5]	500	µg/g
Tin (Sn)	[7440-31-5]	0,05	%
Tin (Sn)	[7440-31-5]	500	µg/g
Titanium (Ti)	[7440-32-6]	0,05	%
Titanium (Ti)	[7440-32-6]	500	µg/g
Vanadium (V)	[7440-62-2]	0,05	%
Vanadium (V)	[7440-62-2]	500	µg/g
Zinc (Zn)	[7440-66-6]	0,05	%
Zinc (Zn)	[7440-66-6]	500	µg/g