

**CONOSTAN - Multi-element standard (custom blend) AM-special in blank oil: 750 ppm**

Art. ID                      CON-150-250-026  
Unit                         200 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	0,075	%		
	Barium (Ba)	[7440-39-3]	750	µg/g		
	Calcium (Ca)	[7440-70-2]	0,075	%		
	Calcium (Ca)	[7440-70-2]	750	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,075	%		
	Magnesium (Mg)	[7439-95-4]	750	µg/g		
	Phosphorus (P)	[7723-14-0]	0,075	%		
	Phosphorus (P)	[7723-14-0]	750	µg/g		
	Zinc (Zn)	[7440-66-6]	0,075	%		
	Zinc (Zn)	[7440-66-6]	750	µg/g		