

CONOSTAN - Multi-element standard (custom blend) AM-special in blank oil +B+Mo: 900 ppm

Art. ID CON-150-250-036
Unit 200 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	0,09	%		
	Barium (Ba)	[7440-39-3]	900	µg/g		
	Boron (B)	[7440-42-8]	0,09	%		
	Boron (B)	[7440-42-8]	900	µg/g		
	Calcium (Ca)	[7440-70-2]	0,09	%		
	Calcium (Ca)	[7440-70-2]	900	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,09	%		
	Magnesium (Mg)	[7439-95-4]	900	µg/g		
	Molybdenum (Mo)	[7439-98-7]	0,09	%		
	Molybdenum (Mo)	[7439-98-7]	900	µg/g		
	Phosphorus (P)	[7723-14-0]	0,09	%		
	Phosphorus (P)	[7723-14-0]	900	µg/g		
	Zinc (Zn)	[7440-66-6]	0,09	%		
	Zinc (Zn)	[7440-66-6]	900	µg/g		