

CONOSTAN - Multi-element standard (custom blend) AM-special in blank oil +K: 500 ppm

Art. ID CON-150-250-037
Unit 200 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	0,05	%		
	Barium (Ba)	[7440-39-3]	500	µg/g		
	Calcium (Ca)	[7440-70-2]	0,05	%		
	Calcium (Ca)	[7440-70-2]	500	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,05	%		
	Magnesium (Mg)	[7439-95-4]	500	µg/g		
	Phosphorus (P)	[7723-14-0]	0,05	%		
	Phosphorus (P)	[7723-14-0]	500	µg/g		
	Potassium (K)	[7440-09-7]	0,05	%		
	Potassium (K)	[7440-09-7]	500	µg/g		
	Zinc (Zn)	[7440-66-6]	0,05	%		
	Zinc (Zn)	[7440-66-6]	500	µg/g		