

**CONOSTAN - Multi-element standard (custom blend) AM-special in blank oil +Sb: 3000 ppm**

Art. ID                      CON-150-250-038  
Unit                         200 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	0,3	%		
	Antimony (Sb)	[7440-36-0]	3000	µg/g		
	Barium (Ba)	[7440-39-3]	0,3	%		
	Barium (Ba)	[7440-39-3]	3000	µg/g		
	Calcium (Ca)	[7440-70-2]	0,3	%		
	Calcium (Ca)	[7440-70-2]	3000	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,3	%		
	Magnesium (Mg)	[7439-95-4]	3000	µg/g		
	Phosphorus (P)	[7723-14-0]	0,3	%		
	Phosphorus (P)	[7723-14-0]	3000	µg/g		
	Zinc (Zn)	[7440-66-6]	0,3	%		
	Zinc (Zn)	[7440-66-6]	3000	µg/g		