

## **Reactor steel, No. 1.6310 - chips**

Art. ID                      ECRM-D 192-1  
Unit                         100 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,1875	%		
	Silicon (Si)	[7440-21-3]	0,219	%		
	Manganese (Mn)	[7439-96-5]	1,377	%		
	Phosphorus (P)	[7723-14-0]	0,0029	%		
	Sulfur (S)	[7704-34-9]	0,001	%		
	Chromium (Cr)	[7440-47-3]	0,0717	%		
	Molybdenum (Mo)	[7439-98-7]	0,482	%		
	Nickel (Ni)	[7440-02-0]	0,755	%		
	Aluminium (Al)	[7429-90-5]	0,0308	%		
	Aluminium (Al)	[7429-90-5]	0,0285	%	Acid sol.	
	Cobalt (Co)	[7440-48-4]	0,055	%		
	Copper (Cu)	[7440-50-8]	0,453	%		
	Nitrogen (N)	[7727-37-9]	0,0118	%		
	Alloy / Werkstoff		1.6310			