

**Kit of ethanolic Solutions 9%, 11.50%, 14% (ISO 17034 \ ISO 17025)**

Art. ID	LM24-OP-2004075
Unit	3x100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

**Description**

Standards of Ethanol in Water / These standards are intended for use as reference solutions for alcohol content / measuring instruments. / They are prepared using organic-free water, and the certified concentrations are based on results obtained by density determination. / Accuracy is  $\pm 0.02\%$