

Kit of ethanolic Solutions 50%, 60%, 70%, 80% (ISO 17034 \ ISO 17025)

Art. ID	LM24-OP-2004078
Unit	4 x 100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

Standards of Ethanol in Water / These standards are intended for use as reference solutions for alcohol content / measuring instruments. / They are prepared using organic-free water, and the certified concentrations are based on results obtained by density determination. / Accuracy is $\pm 0.02\%$