

Stainless steel, 13.7% chromium no. 1.4034 - chips

Art. ID ECRM-D 226-1
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,416	%		
	Silicon (Si)	[7440-21-3]	0,514	%		
	Manganese (Mn)	[7439-96-5]	0,434	%		
	Phosphorus (P)	[7723-14-0]	0,0207	%		
	Sulfur (S)	[7704-34-9]	0,0094	%		
	Chromium (Cr)	[7440-47-3]	1,67	%		
	Molybdenum (Mo)	[7439-98-7]	0,024	%		
	Nickel (Ni)	[7440-02-0]	0,139	%		
	Nitrogen (N)	[7727-37-9]	0,0362	%		
	Vanadium (V)	[7440-62-2]	0,022	%		
	Alloy / Werkstoff		1.4034			