

Austenitic manganese steel, 12,7% manganese, no. 1.3401 - chips

Art. ID ECRM-D 235-1
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,912	%		
	Silicon (Si)	[7440-21-3]	0,094	%		
	Manganese (Mn)	[7439-96-5]	12,73	%		
	Phosphorus (P)	[7723-14-0]	0,045	%		
	Sulfur (S)	[7704-34-9]	0,0072	%		
	Chromium (Cr)	[7440-47-3]	0,354	%		
	Copper (Cu)	[7440-50-8]	0,073	%		
	Molybdenum (Mo)	[7439-98-7]	0,032	%		
	Nitrogen (N)	[7727-37-9]	0,02	%		
	Alloy / Werkstoff		1.3401			