

Starch TS (US Pharmacopoeia)

Art. ID	LM24-OP-2004647
Unit	100 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Test Solution (TS) USP, Shelf-life: 6 months /// The reagent is specially intended for use in testing USP and monographs. This solution was produced according USP specifications using USP reagent quality components and water meeting the USP monograph requirement for Purified Water /// Instruction for use: The solution bottle should be open for the minimum time required to dispense the solution. After use, the bottle should be tightly recapped and stored under normal laboratory conditions