

**Ethanol 90% (v/v) (ISO 17034 \ ISO 17025)**

Art. ID	LM24-OP-2004791
Unit	500 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1170 II E2

**Description**

Standards of Ethanol in Water / These standards are intended for use as reference solutions for alcohol content / measuring instruments. / They are prepared using organic-free water, and the certified concentrations are based on results obtained by density determination. / Accuracy is  $\pm 0.02\%$