

Hydrochloric Acid 0.2M (0.2N) (ISO 17034 \ ISO 17025)

Art. ID	LM24-OP-2004863
Unit	500 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

"Analytical Volumetric Solutions / Titration is the oldest but still the most precise procedure in chemical analysis. / However, such reliable analyses require accurately adjusted volumetric solutions. / Produced and calibrated acc. to the instructions given in ISO Guide 34 and ISO 17025 / Traceability to SI of NIST (SRM) or of Laboratories/ Producers accredited according to ISO/ IEC 17025 and/ or ISO Guide 34 / Certificate of analysis prepared acc. to ISO Guide 31:2000 ""Reference materials - Contents of certificates and labels"" /// Produced under ISO 17034 and ISO 17025."