

High-alloyed steel (toolsteel 1.2344) - disc (37 mm x 25 mm)

Art. ID ECRM-D 271-1-D
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,3698	%		
	Silicon (Si)	[7440-21-3]	0,923	%		
	Manganese (Mn)	[7439-96-5]	0,437	%		
	Phosphorus (P)	[7723-14-0]	0,012	%		
	Sulfur (S)	[7704-34-9]	0,00045	%		
	Chromium (Cr)	[7440-47-3]	5,002	%		
	Molybdenum (Mo)	[7439-98-7]	1,247	%		
	Nickel (Ni)	[7440-02-0]	0,1552	%		
	Aluminium (Al)	[7429-90-5]	0,0234	%		
	Arsenic (As)	[7440-38-2]	0,0057	%		
	Cobalt (Co)	[7440-48-4]	0,0139	%		
	Copper (Cu)	[7440-50-8]	0,1371	%		
	Nitrogen (N)	[7727-37-9]	0,0137	%		
	Tin (Sn)	[7440-31-5]	0,0084	%		
	Titanium (Ti)	[7440-32-6]	0,002	%		
	Vanadium (V)	[7440-62-2]	0,85	%		
	Tungsten (W)	[7440-33-7]	0,0054	%		
	Calcium (Ca)	[7440-70-2]	0,0009	%		
	Oxygen (O)	[7782-44-7]	0,002	%		
	Alloy / Werkstoff		1.2344			