

**Potassium Hydroxide in ethanol 0.1M (0.1N) (ISO 17034 \ ISO 17025)**

|                 |   |
|-----------------|---|
| Art. ID         | LM24-OP-2004983                                     |
| Unit            | 500 mL  |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) |

**Description**

"Analytical Volumetric Solutions / Titration is the oldest but still the most precise procedure in chemical analysis. / However, such reliable analyses require accurately adjusted volumetric solutions. / Produced and calibrated acc. to the instructions given in ISO Guide 34 and ISO 17025 / Traceability to SI of NIST (SRM) or of Laboratories/ Producers accredited according to ISO/ IEC 17025 and/ or ISO Guide 34 / Certificate of analysis prepared acc. to ISO Guide 31:2000 ""Reference materials - Contents of certificates and labels"" /// Produced under ISO 17034 and ISO 17025."