

Potassium Hydroxide in ethanol 0.5M (0.5N) (ISO 17034 \ ISO 17025)

Art. ID LM24-OP-2004993
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2924 II E2

Description

"Analytical Volumetric Solutions / Titration is the oldest but still the most precise procedure in chemical analysis. / However, such reliable analyses require accurately adjusted volumetric solutions. / Produced and calibrated acc. to the instructions given in ISO Guide 34 and ISO 17025 / Traceability to SI of NIST (SRM) or of Laboratories/ Producers accredited according to ISO/ IEC 17025 and/ or ISO Guide 34 / Certificate of analysis prepared acc. to ISO Guide 31:2000 ""Reference materials - Contents of certificates and labels"" /// Produced under ISO 17034 and ISO 17025."

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium Hydroxide	[1310-58-3]				