

## **ICP Multi-Element Standard - Quality Control Standard 2R - 7 components in 5 %HNO<sub>3</sub>**

Art. ID LM24-OP-4010076

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Silver (Ag)	[7440-22-4]	100	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		