

Tuning Standard (XXIV) - 15 components (Standard Equivalent to Merck)

Art. ID LM24-OP-4000006
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Potassium (K) | [7440-09-7] | 500 | mg/L | | |
| | Aluminium (Al) | [7429-90-5] | 50 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 50 | mg/L | | |
| | Barium (Ba) | [7440-39-3] | 50 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 50 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 50 | mg/L | | |
| | Cobalt (Co) | [7440-48-4] | 50 | mg/L | | |
| | Copper (Cu) | [7440-50-8] | 50 | mg/L | | |
| | Lead (Pb) | [7439-92-1] | 50 | mg/L | | |
| | Manganese (Mn) | [7439-96-5] | 50 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 50 | mg/L | | |
| | Nickel (Ni) | [7440-02-0] | 50 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 50 | mg/L | | |
| | Strontium (Sr) | [7440-24-6] | 50 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 50 | mg/L | | |