

Tuning Standard (XXIV) - 15 components (Standard Equivalent to Merck) in HNO₃ 1%

Art. ID LM24-OP-4000007
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]	500	mg/L		
	Aluminium (Al)	[7429-90-5]	50	mg/L		
	Arsenic (As)	[7440-38-2]	50	mg/L		
	Barium (Ba)	[7440-39-3]	50	mg/L		
	Cadmium (Cd)	[7440-43-9]	50	mg/L		
	Chromium (Cr)	[7440-47-3]	50	mg/L		
	Cobalt (Co)	[7440-48-4]	50	mg/L		
	Copper (Cu)	[7440-50-8]	50	mg/L		
	Lead (Pb)	[7439-92-1]	50	mg/L		
	Manganese (Mn)	[7439-96-5]	50	mg/L		
	Molybdenum (Mo)	[7439-98-7]	50	mg/L		
	Nickel (Ni)	[7440-02-0]	50	mg/L		
	Selenium (Se)	[7782-49-2]	50	mg/L		
	Strontium (Sr)	[7440-24-6]	50	mg/L		
	Zinc (Zn)	[7440-66-6]	50	mg/L		