

ICP Single-Element Standard - Interference Check Standard 3 - Sb in 2 %HNO₃

Art. ID LM24-OP-4000104
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Antimony (Sb) | [7440-36-0] | 500 | mg/L | | |