

ICP-MS MIX 8 components in 2% HNO₃

Art. ID LM24-OP-4000118
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	2	mg/L		
	Cadmium (Cd)	[7440-43-9]	2	mg/L		
	Chromium (Cr)	[7440-47-3]	3	mg/L		
	Copper (Cu)	[7440-50-8]	10	mg/L		
	Iron (Fe)	[7439-89-6]	4	mg/L		
	Nickel (Ni)	[7440-02-0]	2	mg/L		
	Lead (Pb)	[7439-92-1]	5	mg/L		
	Zinc (Zn)	[7440-66-6]	10	mg/L		