

## **Ion Chromatography Standard Standard- 3 ions in Water**

Art. ID LM24-OP-4000133  
Unit 500 mL  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphate (PO <sub>4</sub> <sup>3-</sup> )		20	mg/L		
	Chloride (Cl <sup>-</sup> )		0,2	mg/L		
	Sulfate (SO <sub>4</sub> <sup>2-</sup> )		0,2	mg/L		