

**FID Sample - 218 mg/l each of in n-Hexane, 3x1 ml**

Art. ID LM24-OP-4000174  
Unit 3 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Tetradecane	[629-59-4]	218	mg/l		
	n-Pentadecane	[629-62-9]	218	mg/l		
	n-Hexadecane	[544-76-3]	218	mg/l		