

**Tuning Solution for ICP-MS 6 components: 1 µg/L each in HNO3 2% (Standard equivalent to Agilent/Varian)**

Art. ID LM24-OP-4000182  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 8 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium (Li)	[7439-93-2]	1	µg/L		
	Magnesium (Mg)	[7439-95-4]	1	µg/L		
	Yttrium (Y)	[7440-65-5]	1	µg/L		
	Cerium (Ce)	[7440-45-1]	1	µg/L		
	Thallium (Tl)	[7440-28-0]	1	µg/L		
	Cobalt (Co)	[7440-48-4]	1	µg/L		