

PA Tuning Solution 2 - 8 components (Standard equivalent to Agilent)

Art. ID LM24-OP-4000186

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Germanium (Ge)	[7440-56-4]	10	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Palladium (Pd)	[7440-05-3]	10	mg/L		
	Ruthenium (Ru)	[7440-18-8]	10	mg/L		
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Iridium (Ir)	[7439-88-5]	5	mg/L		
	Titanium (Ti)	[7440-32-6]	5	mg/L		