

Internal standard mix for ICP-MS systems - 8 components (Standard equivalent to Agilent)

Art. ID LM24-OP-4000187
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium-6 (LI-6)	[14258-72-1]	10	mg/L		
	Scandium (Sc)	[7440-20-2]	10	mg/L		
	Germanium (Ge)	[7440-56-4]	10	mg/L		
	Rhodium (Rh)	[7440-16-6]	10	mg/L		
	Indium (In)	[7440-74-6]	10	mg/L		
	Terbium (Tb)	[7440-27-9]	10	mg/L		
	Lutetium (Lu)	[7439-94-3]	10	mg/L		
	Bismuth (Bi)	[7440-69-9]	10	mg/L		