

6020 Interference Check Solution A for ICP-MS systems - 12 components in Nitric Acid 5% Hydrofluoric acid tr% (Standard equivalent to Agilent)

Art. ID LM24-OP-4000188
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Chloride (Cl-) | | 20000 | mg/L | | |
| | Calcium (Ca) | [7440-70-2] | 3000 | mg/L | | |
| | Iron (Fe) | [7439-89-6] | 2500 | mg/L | | |
| | Sodium (Na) | [7440-23-5] | 2500 | mg/L | | |
| | Carbon (C) | [7440-44-0] | 2000 | mg/L | | |
| | Aluminium (Al) | [7429-90-5] | 1000 | mg/L | | |
| | Magnesium (Mg) | [7439-95-4] | 1000 | mg/L | | |
| | Phosphorus (P) | [7723-14-0] | 1000 | mg/L | | |
| | Potassium (K) | [7440-09-7] | 1000 | mg/L | | |
| | Sulfur (S) | [7704-34-9] | 1000 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 20 | mg/L | | |
| | Titanium (Ti) | [7440-32-6] | 20 | mg/L | | |