

6020 Interference Check Solution B for ICP/MS - 11 components in HNO3 5% (Standard equivalent to Agilent)

Art. ID LM24-OP-4000190
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Cobalt (Co) | [7440-48-4] | 20 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 20 | mg/L | | |
| | Copper (Cu) | [7440-50-8] | 20 | mg/L | | |
| | Manganese (Mn) | [7439-96-5] | 20 | mg/L | | |
| | Nickel (Ni) | [7440-02-0] | 20 | mg/L | | |
| | Vanadium (V) | [7440-62-2] | 20 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 10 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 10 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 10 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 10 | mg/L | | |
| | Silver (Ag) | [7440-22-4] | 5 | mg/L | | |