

ICP-MS tuning solution- 6 components in HNO3 2% (Standard equivalent to Agilent/ Varian Equivalent)

Art. ID LM24-OP-4000192
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium (Li)	[7439-93-2]	10	µg/mL		
	Magnesium (Mg)	[7439-95-4]	10	µg/mL		
	Yttrium (Y)	[7440-65-5]	10	µg/mL		
	Cerium (Ce)	[7440-45-1]	10	µg/mL		
	Thallium (Tl)	[7440-28-0]	10	µg/mL		
	Cobalt (Co)	[7440-48-4]	10	µg/mL		