

**ICP Multi-Element Standard - EPA 200.7 - Mixed Calibration Standard 4 - 5 components in 5 %HNO<sub>3</sub>**

Art. ID LM24-OP-4000196

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	50	mg/L		
	Zinc (Zn)	[7440-66-6]	50	mg/L		
	Tin (Sn)	[7440-31-5]	40	mg/L		