

ICP-OES wavelength calibration solution - 6 components (Standard equivalent to Agilent)

Art. ID LM24-OP-4000211
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	100	mg/L		
	Indium (In)	[7440-74-6]	100	mg/L		
	Lithium-6 (LI-6)	[14258-72-1]	100	mg/L		
	Scandium (Sc)	[7440-20-2]	100	mg/L		
	Terbium (Tb)	[7440-27-9]	100	mg/L		
	Yttrium (Y)	[7440-65-5]	100	mg/L		