

Calibration mix 1 for AA and ICP-OES - 4 components (Standard equivalent to Agilent)

Art. ID LM24-OP-4000212
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		