

**Calibration mix majors (Standard equivalent to Agilent) 500 mg/L each of in 5% HNO<sub>3</sub>**

Art. ID LM24-OP-4000214

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	500	mg/L		
	Iron (Fe)	[7439-89-6]	500	mg/L		
	Potassium (K)	[7440-09-7]	500	mg/L		
	Magnesium (Mg)	[7439-95-4]	500	mg/L		
	Sodium (Na)	[7440-23-5]	500	mg/L		