

Calibration mix majors (Standard equivalent to Agilent)

Art. ID LM24-OP-4000215
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Calcium (Ca) | [7440-70-2] | 500 | mg/L | | |
| | Iron (Fe) | [7439-89-6] | 500 | mg/L | | |
| | Potassium (K) | [7440-09-7] | 500 | mg/L | | |
| | Magnesium (Mg) | [7439-95-4] | 500 | mg/L | | |
| | Sodium (Na) | [7440-23-5] | 500 | mg/L | | |