

Multi-element calibration standard-3 - 10 components (Standard equivalent to Agilent)

Art. ID LM24-OP-4000244
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Gold (Au)	[7440-57-5]	10	mg/L		
	Hafnium (Hf)	[7440-58-6]	10	mg/L		
	Iridium (Ir)	[7439-88-5]	10	mg/L		
	Palladium (Pd)	[7440-05-3]	10	mg/L		
	Platinum (Pt)	[7440-06-4]	10	mg/L		
	Rhodium (Rh)	[7440-16-6]	10	mg/L		
	Ruthenium (Ru)	[7440-18-8]	10	mg/L		
	Tellurium (Te)	[13494-80-9]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		