

ICP Multi-Element Standard - ICP Calibration Standard - Toxic components - 7 components in 5 %HNO₃

Art. ID LM24-OP-4000249

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		