

**6020 Interference Check Solution A for ICP-MS systems 12 components in Nitric Acid 5% Hydrofluoric acid tr% (Standard equivalent to Agilent/Varian)**

Art. ID LM24-OP-4000258  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chloride (Cl-)		20000	µg/mL		
	Calcium (Ca)	[7440-70-2]	3000	µg/mL		
	Iron (Fe)	[7439-89-6]	2500	µg/mL		
	Sodium (Na)	[7440-23-5]	2500	µg/mL		
	Carbon (C)	[7440-44-0]	2000	µg/mL		
	Aluminium (Al)	[7429-90-5]	1000	µg/mL		
	Magnesium (Mg)	[7439-95-4]	1000	µg/mL		
	Phosphorus (P)	[7723-14-0]	1000	µg/mL		
	Potassium (K)	[7440-09-7]	1000	µg/mL		
	Sulfur (S)	[7704-34-9]	1000	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	20	µg/mL		
	Titanium (Ti)	[7440-32-6]	20	µg/mL		