

**6020 Interference Check Solution B for ICP/MS - 11 components in HNO3 5% (Standard equivalent to Agilent/Varian)**

Art. ID LM24-OP-4000263  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	20	µg/mL		
	Cobalt (Co)	[7440-48-4]	20	µg/mL		
	Copper (Cu)	[7440-50-8]	20	µg/mL		
	Manganese (Mn)	[7439-96-5]	20	µg/mL		
	Nickel (Ni)	[7440-02-0]	20	µg/mL		
	Vanadium (V)	[7440-62-2]	20	µg/mL		
	Arsenic (As)	[7440-38-2]	10	µg/mL		
	Cadmium (Cd)	[7440-43-9]	10	µg/mL		
	Selenium (Se)	[7782-49-2]	10	µg/mL		
	Zinc (Zn)	[7440-66-6]	10	µg/mL		
	Silver (Ag)	[7440-22-4]	5	µg/mL		