

Arsenic (As) 1000 mg/L in diluted HNO₃, Single-Element Standard for AAS (Flame)

Art. ID LM24-OP-4000280
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Arsenic (As) | [7440-38-2] | 1000 | mg/L | | |