

**Tin (Sn) 1000 mg/L in 2-5% HNO<sub>3</sub>/ tr.HF, Single-Element Standard for AAS (Flame)**

Art. ID LM24-OP-4000423  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	1000	mg/L		