

**Beryllium (Be) 10000 mg/L in 2-5% HNO<sub>3</sub>/ tr.HF, Single-Element Concentrate Standard for AAS (Flame)**

Art. ID LM24-OP-4000469  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	10000	mg/L		