

**Tin (Sn) 10000 mg/L in 2-5% HNO<sub>3</sub>/ tr.HF, Single-Element Concentrate Standard for AAS (Flame)**

Art. ID LM24-OP-4000505  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	10000	mg/L		